




SURGICAL SEQUENCE

	PILOT DRILL	STOP DRILL N°1	DRILL N°2	DRILL N°3	DRILL N°4	DRILL N°5	SCREW TAP Ø3,5	SCREW TAP Ø4	SCREW TAP Ø5	COUNTERSINK Ø3,5	COUNTERSINK Ø4	COUNTERSINK Ø5
REF.	A-FPT310	A-FBXxxx <i>(depending on implant length)</i>	A-FCX200	A-FCX300	A-FCX400	A-FCX500	A-TAR306	A-TAR406	A-TAR506	A-ALE300	A-ALE400	A-ALE500
RPM*	1200	1200	1000	800	600	500	15	15	15	500 to 600	500 to 600	500 to 600
Ø3,5	●	●	●	▶●			▶●			▶●		
Ø4	●	●	●	●	▶●			▶●			▶●	
Ø5	●	●	●	●	●	▶●			▶●			▶●
	Trephine the cortical bone with the pilot drill to facilitate the penetration of the first drill (1200 rpm)*.	Use the stop drill n°1 fitted to the length of the implant (1200 rpm)*.	Use drill n°2 to the required length (1000 rpm)*.	For the implants Ø3.5; Ø4 and Ø5: use the drill n°3 to the required length (800 rpm)*.	For the implants Ø4 and Ø5: use the drill n°4 to the required length (600 rpm)*.	For the implants Ø5: use the drill n°5 to the required length (500 rpm)*.	Use the screw tap Ø3.5 for the implants Ø3.5 to the required length (15 rpm)*.	Use the screw tap Ø4 for the implants Ø4 to the required length (15 rpm)*.	Use the screw tap Ø5 for the implants Ø5 to the required length (15 rpm)*.	Use the countersink corresponding to the diameter the zirconia collar. Ream the cortical bone up to the laser marking (500 to 600 rpm)*.		

*The rotation speeds indicated are for information only and depend on the bone quality.

○ Use of the drill depending on the length of the implant to be placed ● For Ø3.5 implants ● For Ø4 implants ● For Ø5 implants



SCANORA AND X-RAY TEMPLATE:

Product code: A-TS600

The Z1-Connect Implant (diameter and length) is selected using the scanora and X-ray template.

Take into account the tip of the drills which is 1mm long while evaluating the available bone height.









SURGICAL KIT:







Product code: A-TCP006

All the instruments needed to place the Z1-Connect Implants are available in the TBR surgical kit.

CONTRA-ANGLE PROTOCOL

	CONTRA-ANGLE SCREWTOOL			MANUAL SCREWTOOL		
REF.	A-MCA325 [long] - A-MCA230 [short]			A-MCC258 [long] - A-MCC163 [short]		
RPM	N/A	15 to 20		Manual		
	 <p>Take the screwdriver for contra-angle and clamp the implant inside its packaging. Maintain the contra-angle facing up while moving the implant to the surgical site.</p>	 <p>Screw the implant in the alveolar ridge until the implant is completely inserted. <i>NB: Index the implant connection using the visual mark on the screwdriver. The laser marking indicates the position of a side of the internal octagon of the implant connection.</i></p>	 <p>Remove the contra-angle vertically as well as the contra-angle screwdriver.</p>	 <p>If ever the implant is not completely screwed in using the contra-angle, finish the insertion with the torque-ratchet wrench [GAN-469-1000203] and its screwdriver.</p>	 <p>Remove the cover screw from its packaging using the screwdriver for torque-ratchet wrench [GAN-469-1000203]. Maintain the screwdriver pointing up while transporting the screw to the surgical site. Seal the implant with the cover screw.</p>	 <p>Suture the gum. Check radiologically that the implant is perfectly positioned in the bone.</p>

TORQUE-RATCHET WRENCH PROTOCOL

	MANUAL SCREWTOOL					
REF.	A-MCC258 [long] - A-MCC163 [short]					
RPM	N/A	Manual				
	 <p>Take the screwdriver and clamp the implant inside its packaging. Maintain the screwdriver pointing up while transporting the implant to the surgical site.</p>	 <p>Begin screwing the implant manually.</p>	 <p>Finish screwing with the torque-ratchet wrench [GAN-469-1000203]. Screw the implant in the alveolar ridge until the implant is completely inserted. <i>NB: Index the implant connection using the visual mark on the screwdriver. The laser marking indicates the position of a side of the internal octagon of the implant connection.</i></p>	 <p>Remove the torque-ratchet wrench [GAN-469-1000203] and pull the screwdriver out vertically.</p>	 <p>Remove the cover screw from its packaging. Maintain the screwdriver pointing up while transporting the screw to the surgical site. Seal the implant with the cover screw.</p>	 <p>Suture the gum. Check radiologically that the implant is perfectly positioned in the bone.</p>